



# ToxCO<sup>®</sup>

*Saving lives, one breath at a time.*

## Features and Benefits




For use with adults



To alert staff of high atmospheric CO



For use with unconscious patients



Tag Readings



Maintenance Reminders



Easy to use touchscreen interface

## Ideal for:

- General practice
- Ambulance services
- Firefighters
- First responders
- Health visitors



## Technical Specification

Breath test Concentration range (CO)	0-50%COHb/0-500ppm
Face mask test Concentration range (CO)	0-28%COHb/0-200ppm
Ambient Test Concentration range (CO)	0-500ppm
Display	Full colour touchscreen
Detection principle	Electrochemical sensor
Repeatability	≤±5% difference on consecutive readings
Accuracy	≤±3ppm/10% - whichever is greater*
Power	3 x AA (LR6 or equivalent) - up to 1000 minutes 1 x CR2032 lithium coin cell
T <sub>90</sub> response time	<30 seconds
Operating temperature	0-45°C
Storage/ transport temperature	0-50°C
Operating/storage/transport pressure	Atmospheric ±10% (912-1114mbar)
Operating humidity	15-90% non-condensing
Storage/transport humidity	5-95%
Sensor operating life	2 years
Sensor sensitivity	1ppm
Sensor drift	<5% per annum
Dimensions	Approx. 37 x 77 x 140mm
Weight	Approx. 215g (including batteries)
Materials	Case: polycarbonate/ABS blend with Antimicrobial additive D-piece™: polypropylene SteriBreath™ Eco: paper OneBreath™: polypropylene
H <sub>2</sub> cross interference	≤6%

\*Readings of >500ppm at temperatures between 0-14°C can decrease accuracy to ≤±3ppm/15%.

Visit [www.bedfont.com/resources](http://www.bedfont.com/resources) to view this document in other languages.



Bedfont® Scientific Ltd.  
Station Road, Harrietsham, Maidstone,  
Kent, ME17 1JA England  
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860  
Email: [ask@bedfont.com](mailto:ask@bedfont.com) Web: [www.bedfont.com](http://www.bedfont.com)



Emergo Europe B.V.  
Westervoortsedijk 60  
6827 AT Arnhem  
The Netherlands.

© Bedfont® Scientific Limited 2024

Issue 13 - September 2024, Part No: MKT376  
Bedfont® Scientific Limited reserves the right to change or update his literature without prior notice.  
Registered in: England and Wales. Registered No: 1289798



MID 502905

