


CO testing with the iCOquit[®] Smokerlyzer[®] couldn't be easier:

1

Start a breath test in the App



2

Inhale and hold breath



3

Exhale into the iCOquit[®]



4

Read your result



For more information
on the iCOquit[®]
scan here
Or visit:



iCOquit[®] Smokerlyzer[®]

For use with the iCOquit[®] app
and the iCOquit[®] Mums to Be app



Bedfont[®] Scientific Ltd.
Station Road, Harrietsham, Maidstone,
Kent, ME17 1JA England
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860
Email: ask@bedfont.com Web: www.bedfont.com



Emergo Europe B.V.
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands.

© Bedfont[®] Scientific Limited 2024

Issue 3 - September 2024, Part No: MKT561

Bedfont[®] Scientific Limited reserves the right to change or update this literature without prior notice. Registered in: England and Wales. Registered No: 1289798

